


Intensified Pumps for Intensified Processes

RSC Symposium
June 6th 2013, Munich

Dr. Carsten Damerou

HNP Mikrosysteme GmbH
Bleicherufer 25
D-19053 Schwerin
Telefon +49 385 52190-347
Telefax +49 385 52190-333
E-mail carsten.damerou@hnp-mikrosysteme.de
Internet www.hnp-mikrosysteme.de



- **Introduction**
- mzs-pumps
- **Characteristic properties**
- **R&D, materials**
- **Example**

Pump classification



Liquid pumps

Positive displacement pumps

Rotary pumps

- External gear pump
- Internal gear pump
- Annular gear pump
- Lobe pump
- Vane pump
- Rotary piston pump
- Impeller pump
- Spiral pump
- Screw pump
- Eccentric screw pump
- Peristaltic pump

Reciprocating pumps

- Diaphragm pump
- Piston pump
- Syringe Pump
- Reciprocating plunger pump
- Plunger pump

Rotodynamic pumps

Rotary pumps

- Axial-flow pump
- Centrifugal pump
- Side channel pump
- Peripheral pump

Reciprocating pumps

Gas pumps

- compressor
- vacuum pump
- blower

rotary

- Axial blower
- Vane pump
- Screw compressor
- Roots compressor
- Scroll compressor

reciprocating

- Diaphragm pump
- Piston compressor

miscellaneous

- Eductor-jet pump
- Cryogenic pump
- Sorption pump

Water jet pump

RSC_SympChemspec_2013-06-06

A common enquiry!

Flow rate range


1 μ l/h – 1 l/min

Viscosity range

0.2 – 1,000,000 mPas

Liquids

liquid butane ...
nitric acid



Differential pressure

Vacuum, 0 – 200 bar

Dosage

2 μ l-droplets,
100 per minute

Dosage precision

< 1 %

Flow

pulseless

- Introduction
- **mzr-pumps**
- Characteristic properties
- R&D, materials
- Example

Micro annular gear technology

Milestones R & D



1996

hydraulic motor diameter 10 mm



1997

micro annular gear pump
outer diameter 2,5 mm



1998

high performance series
mzr-2900, mzr-4600



1999

integrated mzr-pump



2000

mzr-pump with universal actuator



2001

low pressure pump mzr-2921

- cycloidal gear technology
- tungsten carbide, ceramics
- wear-resistant
- long service life
- high contour precision $\pm 1 \mu\text{m}$
- wire EDM, grinding, lapping

Gear type

Diameter of outer rotor
Displacement volume

mzr-2500
Ø 3,4 mm
1,5 μl

mzr-6300
Ø 9,0 mm
24 μl

mzr-2900
Ø 3,4 mm
3 μl


mzr-7200
Ø 9,0 mm
48 μl

mzr-4600
Ø 5,4 mm
12 μl

mzr-11500
Ø 14,0 mm
192 μl


Design and basic principle

Micro annular gear pump




Basic principle


Micro annular gear pump




0°, 309°
0°, 360°




39°
45°




77°
90°




116°
135°




154°
180°



193°
225°




231°
270°



270°
315°

Clearance Space

Comparison of dimensions



human hair





fog droplet



red blood corpuscle


rotor dimensions in natural scale



clearance space
micro annular gear pumps

- Introduction
- mzs-pumps
- **Characteristic properties**
- R&D, materials
- **Example**


Comparison of pulsation



Pulsation of flow:

$$\delta_Q = \frac{Q_{\max} - Q_{\min}}{Q}$$


Diaphragm Pump



$$\delta_Q = \frac{3,14 - 0}{1} = 314 \%$$


Factor 50

Micro annular gear pump



mzr-2921
6/7 indenting

$$\delta_Q = \frac{1,02 - 0,96}{1} = 6,3 \%$$



mzr-2521
10/11 indenting

$$= 1,5 \%$$

Comparison of pulsation – experimental




Tabelle 5.2: Übersicht der Messungen der Flussratengenauigkeit von verschiedenen Pumpenmodellen.

	MZR	Cavro	Cavro (Pendel- betrieb)	Harvard	Vici	Schlauch
mittlere Flussrate [$\mu\text{L}/\text{min}$]	110.2	100.8	101.8	101.0	106.3	98.3
Standard- abweichung [%]	2.5	1.5	14.1	2.7	1.7	33.6
F_{min} [$\mu\text{L}/\text{min}$]	102.5	97.8	41.7	93.5	99.3	35.4
ΔF_{min} [%]	7.0	3.0	58.8	7.5	6.6	64.0
F_{max} [$\mu\text{L}/\text{min}$]	115.6	102.9	174.0	107.9	113.1	141.8
ΔF_{max} [%]	4.9	2.1	72.1	6.8	6.4	44.2

Abb. 5.2: Mit einem Flusssensor wurde die Pulsation verschiedener Pumpentypen bei einer Flussrate von $\sim 100 \mu\text{L}/\text{min}$ mit 1.5625 Hz aufgezeichnet. Die höchste Pulsation zeigte erwartungsgemäß die Schlauchpumpe. Einen nahezu pulsationsfreien Fluss generierte die Kolbenpumpe der Firma Vici. Die hohen Ausschläge der Spritzenpumpe von Cavro sind auf den Pendelbetrieb zurückzuführen.


Source: Fabel, Susanne; PhD thesis, Munich 2007, p. 97

Dosage Precision

Dispensing with mzs-pumps



CV values achievable at 1 mPas



Definition

The coefficient of variation (CV) expressed in % is defined as a standard deviation from the average value.

$$CV [\%] = \frac{s \cdot 100}{\bar{V}}$$

CV Coefficient of variation [%]

s Standard deviation [µl]

V Average volume value [µl]

Remarks

- The CV values have been established for water (1 mPas) by using an analytical scale with 0,01 mg resolution
- Characteristics curves show values obtained at optimal conditions for each individual pump size
- Precision improves with increasing viscosity and decreasing differential pressure


Requirements for optimal precision

- Degassed liquid
- Absence of air bubbles in the system
- Constant pressure conditions
- Return valves preventing free flow
- Well adapted dosing nozzle

- Introduction
- mzs-pumps
- Characteristic properties
- **R&D, materials**
- Example

Testing of novel materials

Long term wear test



Test condition
- difficult liquid
- delta p: 20 bar
- 20 °C
- 3000 rpm

- Introduction
- mzs-pumps
- Characteristic properties
- R&D, materials
- **Example**

mzr® micro annular gear pumps

Product summary



Ex-Pumps

- for pump heads of high performance series and hermetic inert series
- Ex-certification ATEX, EU directive 94/9/EEC
- CE Ex II 2 G c T4 X, CE Ex II 2 G c T5 X

High performance series

- industrial equipment
- tungsten carbide Ni-based, stainless steel 316L seals: PTFE, FKM, optional: EPDM, FFKM
- differential pressure range 40 bar (max. 150 bar)
- wide viscosity range (0.3-1,000,000 mPas)
- DC-servomotor with integrated controller
- modular system: fluidic seal module, heat insulation module, electrical heating, double shell heating and cooling module, reduction gear

Hermetic inert series

- chemically inert materials
- Al₂O₃, ZrO₂ ceramics, alloy C22, SSiC, Kalrez®
- hermetic magnetic coupling
- DC-servomotor with integrated controller



Modular Series

- chemically inert, compact dimensions
- configurable materials: ceramics, alloy C276/C22, optional stainless steel 316L, optional PEEK™;
- seals: PTFE, FKM, optional: EPDM, FFKM.
- DC-motor with graphite brushes

Low pressure series

- compact dimensions
- low pressure range (0-5 bar)
- tungsten carbide Ni-based, stainless steel 316L seals: PTFE, FKM, optional: EPDM, FFKM
- low viscosity liquids (0-100 mPas)
- DC-motor with graphite brushes

Customized pumps

- specific liquids
- specific motor
- specific design
- ...

Micro annular gear pumps

Hermetic inert series



Note:
Optional double shell heating
module for mzr-1155x available

- **High chemical resistance**
against oxidizing and reducing media, acids and bases
- **Hermetically sealed**
magnetic coupling (NdFeB)
- **Long service life**
wear-resistant, ultra-hard materials
- **Compact, chemically inert pump head**
alloy C22, SiC, Al₂O₃ and ZrO₂ ceramics, option: tantalum
- **Precision motor and user-friendly control**
dynamic DC-servomotor with encoder and microcontroller, RS-232 or CAN-Bus, analog input 0-10 V, I/O
- **Precise dosage – low pulsation**
rotary micro annular gear technology, no valves



		mzr-6355	mzr-7255	mzr-1155x
Displacement volume	µl	24	48	192
Smallest dosage volume	µl	15	30	100
Flow rate	ml/min	0.024 – 144	0.048 – 288	0.192 – 1152
Max. operation pressure *	bar	80	80	60 (max. 190)
Differential pressure range	(water, 1 mPas)	0 – 15	0 – 20	30
	(oil, 16 mPas)	0 – 40	40	0 – 60
Viscosity range	mPas	0.3 – 1000		
Speed range	rpm	1 – 6000		

* operation pressure = inlet pressure + differential pressure


RSC_SympChemspec_2013-06-06

Micro annular gear pump mZR®-7255


Hermetic inert series



pump head twisted 90° total length: 146 mm



Note: All cross sections and views are schematic!



RSC_SympChemspec_2013-06-06

Thank you for your attention



HNP Mikrosysteme GmbH
Bleicherufer 25
D-19053 Schwerin
Telefon +49 385 52190-300
Telefax +49 385 52190-333
E-mail info@hnp-mikrosysteme.de
Internet www.hnp-mikrosysteme.de